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
Preparation of a 200 fg/ μ L Chloramphenicol Test Solution

1. Prepare a 0.02% acetic acid solution in water:
 - a. Add 200 mL of water to a 250 mL graduated cylinder.
 - b. Add 40 μ L of acetic acid to the same graduated cylinder.
 - c. Mix the contents thoroughly.
2. Prepare a 1 μ g/ μ L chloramphenicol solution:
 - a. Add 10 mg of chloramphenicol to a 10 mL volumetric flask.
 - b. Fill the flask to volume with methanol.
 - c. Mix the contents thoroughly.
3. Prepare a 1 ng/ μ L chloramphenicol solution:
 - a. Use a funnel to add the 0.02% acetic acid solution in water to within one centimeter of the graduation mark of a 100ml volumetric flask.
 - b. Use a 500- μ L syringe to transfer 100 μ L of the 1 μ g/ μ L solution.
 - c. Then, use a pipette to carefully add dropwise the solvent mixture until the meniscus is level with the graduation mark (Figure 1, other side).
 - d. Mix the contents thoroughly.
4. Prepare a 200 fg/ μ L chloramphenicol test solution:
 - a. Use a funnel to add the 0.02% acetic acid solution in water to within one centimeter of the graduation mark of a 100ml volumetric flask.
 - b. Use a 25- μ L syringe to transfer 20 μ L of the 1 ng/ μ L chloramphenicol solution.
 - c. Then, use a pipette to carefully add dropwise the solvent mixture until the meniscus is level with the graduation mark (Figure 1, other side).
 - d. Mix the contents thoroughly.
5. Label the flask:
 - a. Name: TSQ Installation 200 fg/ μ L Chloramphenicol Test Solution
 - b. Mix Date: DDMMYYYY (for example, 26JUN2017)
 - c. Expiration Date: 1 month after the mix date

Preparation of a 200 fg/ μ L Reserpine Test Solution

1. Prepare a 50:50 (v/v) methanol/water mixture containing 0.1% acetic acid:
 - a. Add 100 mL of water and 100 mL methanol to a 250 mL graduated cylinder.
 - b. Open the ampoule containing 1 mL acetic acid.
 - c. Use a 500 μ L syringe to transfer 200 μ L of acetic acid to the graduated cylinder.
 - d. Mix the contents thoroughly.
2. Prepare a 200 fg/ μ L reserpine test solution:
 - a. Use a funnel to add the 50:50 (v/v) methanol/water mixture containing 0.1% acetic acid to within one centimeter of the graduation mark of a 100ml volumetric flask.
 - b. Open the vial containing the 100pg/ μ L reserpine standard.
 - c. Use a 500- μ L syringe to transfer 200 μ L of the reserpine standard.
 - d. Then, use a pipette to carefully add dropwise the solvent mixture until the meniscus is level with the graduation mark (Figure 1).
 - e. Mix the contents thoroughly.
3. Label the flask:
 - a. Name: TSQ Installation 200 fg/ μ L Reserpine Test Solution
 - b. Mix Date: DDMMYYYY (for example, 26JUN2017)
 - c. Expiration Date: 1 month after the mix date

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The screenshot shows the instrument control software interface with several key components labeled in Chinese:

- 开扫描** (Start Scan): Points to the large blue play button.
- 仪器状态设置为standby** (Instrument status set to standby): Points to the pause button.
- 仪器状态设置为Off** (Instrument status set to Off): Points to the stop button.
- 正负扫描模式切换** (Positive/Negative scan mode switch): Points to the 'Positive' mode selector.
- 质谱图采集模式切换** (Mass spectrum acquisition mode switch): Points to the 'Profile' mode selector.
- 废液阀切换** (Waste valve switch): Points to the 'Valve 1-6' selector.
- 谱图平均** (Spectrum averaging): Points to the 'Avg. (2) ON' selector.
- 更改流速** (Change flow rate): Points to the 'Syringe ON' dropdown menu.
- 保存路径** (Save path): Points to the file path input field in the 'Record' dialog.
- 谱图保存名称** (Spectrum save name): Points to the sample name input field in the 'Record' dialog.

Additional visible UI elements include a 'Record' button, a 'View' button, and a 'Scans to Average' dialog box showing the value '2'. The 'Syringe' settings show a 'Flow Rate (µL/min)' of 5 and a 'Volume (µL)' of 500.

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
Signal vs Scan Number

Plot Chromatogram

- Spray Stability
- TIC
- Base Peak
- User Defined m/z

Scan	Mass
1	622

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

System Calibration Report

Date & Time Tuesday, October 10, 2023 12:13:56 PM
Instrument Model TSQ Altis
Instrument Serial TSQ-A-11049
Software Version 3.4.3279

Name	Polarity	Result	Comment
System Calibration	Positive	Successful	

Source Parameters	
Source Type	H-ESI
Spray Voltage (V)	3500
Sheath Gas (Arb)	3
Aux Gas (Arb)	0
Sweep Gas (Arb)	0
Ion Transfer Tube Temp (°C)	325
Vaporizer Temp (°C)	0


System Calibration Report

Date & Time Tuesday, October 10, 2023 11:47:01 AM
Instrument Model TSQ Altis
Instrument Serial TSQ-A-11049
Software Version 3.4.3279

Name	Polarity	Result	Comment
System Calibration	Negative	Successful	

Source Parameters	
Source Type	H-ESI
Spray Voltage (V)	3000
Sheath Gas (Arb)	4
Aux Gas (Arb)	0
Sweep Gas (Arb)	0
Ion Transfer Tube Temp (°C)	325
Vaporizer Temp (°C)	0

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HESI+: Set Detector to Reserpine Test Settings for HESI


HESI-: Set Detector to Chloramphenicol Test Settings for HESI

APCI+: Set Detector to Reserpine Test Settings for APCI

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

Status: PASSED

SNR Limit: 500000

File: HESI positive sensitivity at 0.2 PW.raw

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



Status: PASSED

SNR Limit: 500000

File: HESI negative sensitivity at 0.7 PW.raw

Instrument Serial Number: TSQ-A-11049

Acquisition Date: 10/10/2023


Noise Range RT: 4.52 min - 4.62 min Average: 4.866E-003 Standard Deviation: 1.710E-002
 Peak at RT: 2.32 min Intensity: 11169.80 SNR: 653246.88 Area: 28584.99
 Peak at RT: 4.12 min Intensity: 10899.18 SNR: 637419.94 Area: 28357.25
 Peak at RT: 0.99 min Intensity: 10616.99 SNR: 620916.73 Area: 27006.37

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

Red boxes highlight the flow rate (0.350 [ml/min]), pressure (0 [bar]), and eluent percentages (70.0 [%], 30.0 [%]). A red box around the eluent selection dropdowns is labeled "确认流动相" (Confirm mobile phase).

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Inject control panel with the following settings and buttons:

- Pos: R:E8
- Vol: 5.00 [μl]
- Buttons: Inject, Stop Inject, Cleanup

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
thermo scientific

Signal to Noise Test

Status: PASSED
 SNR Limit: 100000
 File: APCI positive sensitivity at 0.7 PW.raw
 Instrument Serial Number: TSQ-A-11049
 Acquisition Date: 10/10/2023


Noise Range RT: 3.24 min - 3.34 min Average: 2.503E-004 Standard Deviation: 1.371E-003
 Peak at RT: 2.20 min Intensity: 630.97 SNR: 460243.50 Area: 1936.40
 Peak at RT: 3.50 min Intensity: 623.05 SNR: 454470.20 Area: 1968.71
 Peak at RT: 1.13 min Intensity: 597.13 SNR: 435560.50 Area: 1917.64

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The screenshot displays the 'Method Editor' window for a Thermo Xcalibur instrument. The interface includes a 'Method Timeline' at the top, a 'Global Parameters' section with a timeline, and an 'Ion Source Properties' panel on the right. A red box highlights the 'Ion Source Properties' panel, and a red arrow points from the 'Ion Source' icon in the Global Parameters section to this panel.

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Method Editor

Global Parameters Scan Parameters Summary

Method Timeline Experiment ACTIONS Settings

Method Duration (min) 60

Method Timeline: SRM

Experiment # 1 CLEAR

Scans

- SRM
- Full Scan Q1
- Full Scan Q3
- Product Ion Scan
- Precursor Ion Scan
- Neutral Loss Scan
- SIM Q1
- SIM Q3
- QED

Compound	Retention Time (min)	RT Window (min)	Polarity	Precursor (m/z)	Product (m/z)	Collision Energy (V)	Min Dwell Time (ms)
Reserpine	1	1	Positive	609.28	195.1	38	199.271

SRM Properties

- Polarity: Defined in Table
- Use Cycle Time:
- Cycle Time (sec): 0.2
- Use Calibrated RF Lens:
- Q1 Resolution (FWHM): 0.4
- Q3 Resolution (FWHM): 0.7
- CID Gas (mTorr): 2
- Source Fragmentation: 0
- Chromatographic Peak Width (sec): 6
- Use Chromatographic Filter:
- Use Retention Time Reference:
- Display Retention Time:
- Use Quan Ion:
- Show Visualization:

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
The screenshot shows the Yingsheng TraceFinder Clinical LC software interface. The 'File' menu is open, and 'New Main Method' (主方法...) is highlighted. A red arrow points to the 'Create Main Method' dialog box. The dialog box has a title bar '创建主方法' and a close button. The main text reads: '创建主方法有多种方法,选择您要使用的技术。' Below this are several radio button options:

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
At the bottom of the dialog box are '确定' (OK) and '取消' (Cancel) buttons. In the bottom-left corner of the software interface, the 'Method Development' (方法开发) tab is highlighted with a red box.

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

The screenshot displays the Yingsheng TraceFinder Clinical LC software interface. The left sidebar shows the '方法视图' (Method View) menu with '化合物' (Compound) highlighted. The main window shows the '手动峰创建者' (Manual Peak Creator) section. A chromatogram plot is visible with a peak at RT: 1.601. A context menu is open over the peak, with the option '将此峰作为新化合物添加' (Add this peak as a new compound) highlighted. The bottom status bar indicates '没有采集进展' (No acquisition progress).

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The screenshot shows the TraceFinder software interface. A red arrow points from the 'Save' icon in the top toolbar to the 'Save Method' dialog box. The dialog box contains a table of existing methods and a text field for the new method name.

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